Seronorm™
Seronorm™ Lipid
Pathonorm™ L
Pathonorm™ H

The lot-specific document is available at the SERO website

EN

Product description

Seronorm™/Pathonorm™ L/Pathonorm™ H/ Seronorm™ Lipid is a lyophilized animal-based control serum produced from blood collected from veterinary controlled healthy animals at Norwegian facilities which are government approved according to EU regulations. However, since no method can completely exclude the presence of infectious agents, this material should be handled as an ordinary patient sample.

Safety Data Sheet (SDS) is available at our website www.sero.no.

Use of this product as a trueness control is only recommended together with the lot-specific document which includes the analytical data and stability information. Please visit www.sero.no/theserolibrary to download the lot-specific document. The 8-digit number printed on the label at the top of the first page of this document is the code required to access the lot-specific document on the website. If you do not have internet access or the label is missing, please contact your local distributor or SERO directly. SERO's contact information can be found on the last page of this document.

Intended use

Seronorm™/Pathonorm™ L/Pathonorm™ H/Seronorm™ Lipid is intended for use as an assayed quality control material to monitor the precision and trueness of laboratory measurement procedures.

Instructions for use

To secure uniformity, ensure that each vial is similarly treated with respect to preparation for use. Always mix contents before use.

1. Gently tap the vial to ensure the control material is at the bottom of the vial.
2. Remove the screw cap.
3. Lift the rubber stopper without removing it totally. This allows air to enter the vial.
4. Carefully remove the rubber stopper.
5. Add exactly the volume of sterile deionized water, as stated on the vial label.
6. Replace the stopper in the vial and ensure the vial contents are completely dissolved by gentle continuous rolling for approximately 30 minutes. Avoid the formation of foam. DO NOT SHAKE.
7. If the material is to be used for analysis of trace elements, transfer the reconstituted material to a screw cap plastic tube. This is important to prevent contamination from the vial and the rubber stopper.
8. Remove the amount of material required for testing.
9. Promptly replace the stopper in the vial and store the remaining control material at 2-8 °C.
10. Any material removed from the vial but not used should be discarded. DO NOT RETURN control material to the original vial once removed.
11. Establish target values and analytical ranges for the control material based on the measurement procedures in the lab.

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