Seronorm™

Trace Elements Serum

Trace Elements Urine

Trace Elements Whole Blood

The lot-specific document is available at the SERO website

EN

Product description

Seronorm™ Trace Elements Serum, Seronorm™ Trace Elements Urine and Seronorm™ Trace Elements Whole Blood are lyophilized human-based control materials produced from serum, urine and whole blood, respectively, collected from voluntary donors. Each unit is separately checked using CE approved tests to be negative for HBs antigen, HIV p24-antigen and HIV I, II and HCV antibodies. However, since no method can completely exclude the presence of infectious agents, this material should be handled as an ordinary patient sample.

Safety Data Sheet (SDS) is available at our website www.sero.no.

Use of this product as a trueness control is only recommended together with the lot-specific document which includes the analytical data and stability information. Please visit www.sero.no/theserolibrary to download the lot-specific document. The 8-digit number printed on the label at the top of the first page of this document is the code required to access the lot-specific document on the website. If you do not have internet access or the label is missing, please contact your local distributor or SERO directly. SERO's contact information can be found on the last page of this document.

Intended use

Seronorm™ Trace Elements Serum/Seronorm™ Trace Elements Urine/Seronorm™ Trace Elements Whole Blood is intended for use as an assayed quality control material to monitor the precision and trueness of laboratory measurement procedures.

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Instructions for use

To secure uniformity, ensure that each vial is similarly treated with respect to preparation for use. Always mix contents before use.

1. Gently tap the vial to ensure the control material is at the bottom of the vial.
2. Remove the screw cap.
3. Lift the rubber stopper without removing it totally. This allows air to enter the vial.
4. Carefully remove the rubber stopper.
5. Add exactly the volume of sterile deionized water, as stated on the vial label.
6. Replace the stopper in the vial and ensure the vial contents are completely dissolved by gentle continuous rolling for approximately 30 minutes. Avoid the formation of foam. DO NOT SHAKE.
7. Transfer the reconstituted material to a screw cap plastic tube. This is important to prevent contamination from the vial and the rubber stopper.
8. Remove the amount of material required for testing.
9. Promptly replace the stopper in the plastic tube and store the remaining control material at 2-8 °C.
10. Any material removed from the plastic tube but not used should be discarded. DO NOT RETURN control material to the plastic tube once removed.
11. Establish target values and analytical ranges for the control material based on the measurement procedures in the lab.

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