EN

Product description

Seronorm™ Immunoassay Lyo/Seronorm™ Immunoprotein Lyo is a lyophilized human-based control serum produced from blood collected from voluntary blood donors. Each unit is separately checked using CE approved tests to be negative for HBs antigen, HIV p24-antigen and HIV 1, II and HCV antibodies. However, since no method can completely exclude the presence of infectious agents, this material should be handled as an ordinary patient sample.

Safety Data Sheet (SDS) is available at our website www.sero.no.

Use of this product as a trueness control is only recommended together with the lot-specific document which includes the analytical data and stability information. Please visit www.sero.no/theserolibrary to download the lot-specific document. The 8-digit number printed on the label at the top of the first page of this document is the code required to access the lot-specific document on the website. If you do not have internet access or the label is missing, please contact your local distributor or SERO directly. SERO's contact information can be found on the last page of this document.

Intended use

Seronorm™ Immunoassay Lyo/Seronorm™ Immunoprotein Lyo is intended for use as an assayed quality control material to monitor the precision and trueness of laboratory measurement procedures.

Instructions for use

To secure uniformity, ensure that each vial is similarly treated with respect to preparation for use. Always mix contents before use.

1. Gently tap the vial to ensure the control material is at the bottom of the vial.
2. Remove the screw cap.
3. Lift the rubber stopper without removing it totally. This allows air to enter the vial.
4. Carefully remove the rubber stopper.
5. Add exactly the volume of sterile deionized water, as stated on the vial label.
6. Replace the stopper in the vial and place the vial on a roller mixer for approximately 30 minutes until the contents have dissolved completely. Avoid the formation of foam. DO NOT SHAKE.
7. Remove the amount of material required for testing and promptly store the original vial at 2-8 °C. For optimal analyte stability, analyze the material immediately. Discard any remnants of the sample after each use.
8. Prior to subsequent use, remove the vial from 2-8 °C and mix on a roller mixer for a few minutes before sampling. Replace the stopper and promptly store the original vial at 2-8 °C. For optimal analyte stability, analyze the material immediately. Discard any remnants of the sample after each use.
9. DO NOT RETURN control material to the original vial once removed.
10. Establish target values and analytical ranges for the control material based on the measurement procedures in the lab.
11. Refer to the lot-specific product documentation for specific analyte and stability claims.

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