The lot-specific document is available at the SERO website

EN

Product description

Autonorm™ Clin Chem Liq/Autonorm™ Pharmaca Liq is a liquid human-based control serum produced from blood collected from voluntary blood donors. Each unit is separately checked using CE approved tests to be negative for HBs antigen, HIV p24-antigen and HIV I, II and HCV antibodies. However, since no method can completely exclude the presence of infectious agents, this material should be handled as an ordinary patient sample.

Safety Data Sheet (SDS) is available at our website www.sero.no.

The lot-specific document which includes the analytical data and stability information is available online. Please visit www.sero.no/theserolibrary to download the lot-specific document. The 8-digit number printed on the label at the top of the first page of this document is the code required to access the lot-specific document on the website. If you do not have internet access or the label is missing, please contact your local distributor or SERO directly. SERO's contact information can be found on the last page of this document.

Intended use

Autonorm™ Clin Chem Liq/Autonorm™ Pharmaca Liq is intended for in vitro diagnostic use as an assayed quality control material to verify precision of laboratory measurement procedures.

Instructions for use

To secure uniformity, ensure that each vial is similarly treated with respect to preparation for use. Always mix contents thoroughly before use.

- Place the vial on a roller mixer for approximately 30 minutes until the contents have thawed completely.
- Remove the amount of material required for testing and promptly store the original vial at 2-8 °C. For optimal analyte stability, analyze the material immediately. Discard any remnants of the sample after each use. **Bicarbonate in Autonorm™ Clin Chem Liq is extremely volatile. For analysis of Bicarbonate, please follow further instructions described in the lot-specific product documentation.**
- Prior to subsequent use, remove the vial from 2-8 °C and mix on a roller mixer for a few minutes before sampling. Replace the stopper and promptly store the original vial at 2-8 °C. For optimal analyte stability, analyze the material immediately. Discard any remnants of the sample after each use.
- DO NOT RETURN control material to the original vial once removed.
- If the material is to be used for analysis of trace elements, transfer the material to a screw-cap plastic tube. This is important to prevent contamination from the vial and the rubber stopper.
- Establish target values and analytical ranges for the control material based on the measurement procedures in the lab.
- Refer to the lot-specific product documentation for specific analyte and stability claims.

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