The lot-specific document is available at the SERO website

## EN

### Product description

Seronorm™ Urine is a liquid human-based control produced from urine collected from voluntary donors. Each donor is separately checked using CE approved tests to be negative for HBs antigen, HIV p24-antigen and HIV I, II and HCV antibodies. However, since no method can completely exclude the presence of infectious agents, this material should be handled as an ordinary patient sample.

Safety Data Sheet (SDS) is available on request and at our website [www.sero.no](http://www.sero.no).

### Intended use

Seronorm™ Urine is intended for in vitro diagnostic use as an assayed quality control material to verify repeatability and reproducibility of laboratory measurement procedures.

### Instructions for use

1. Place the vial on a roller mixer for approximately 10 minutes.
2. Remove the amount of material required for testing and promptly store the original vial at 2-8 °C. For optimal analyte stability, analyze the material immediately. Discard any remnants of the sample after each use.
3. Prior to subsequent use, remove the vial from 2-8 °C and mix on a roller mixer for a few minutes before sampling. Replace the stopper and promptly store the original vial at 2-8 °C. For optimal analyte stability, analyze the material immediately. Discard any remnants of the sample after each use.
4. DO NOT RETURN control material to the original vial once removed.
5. Establish target values and analytical ranges for the control material based on the measurement procedures in the lab.

### Precautions and warnings

This product contains components classified as follows in accordance with the Regulation (EC) No 1272/2008:

Isothiazolinone chloride

 H317: May cause an allergic skin reaction

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